

ARHGF Polyclonal Antibody

Cat No: HR1AP11279

For research use only

Overview

Product Name	ARHGF Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	ARHGEF15 KIAA0915
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>Rho guanine nucleotide exchange factor 15(ARHGEF15) Homo sapiens Rho GTPases play a fundamental role in numerous cellular processes that are initiated by extracellular stimuli that work through G protein-coupled receptors. This gene encodes a protein that functions as a specific guanine nucleotide exchange factor for RhoA. It also interacts with ephrin A4 in vascular smooth muscle cells. Two alternatively spliced transcripts variants that encode the same protein have been found for this gene. [provided by RefSeq, Aug 2010],</p>