

RHG06 Polyclonal Antibody

Cat No: HR1AP11277

For research use only

Overview

Product Name	RHG06 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	107kDa
GeneID?Human?	ARHGAP6 RHOGAP6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>Rho GTPase activating protein 6(ARHGAP6) Homo sapiens This gene encodes a member of the rhoGAP family of proteins which play a role in the regulation of actin polymerization at the plasma membrane during several cellular processes. This protein is thought to have two independent functions, one as a GTPase-activating protein with specificity for RhoA, and another as a cytoskeletal protein that promotes actin remodeling. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008].</p>