

RHG05 Polyclonal Antibody

Cat No: HR1AP11276

For research use only

Overview

Product Name	RHG05 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	165kDa
GeneID?Human?	ARHGAP5 RHOGAP5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Rho GTPase activating protein 5(ARHGAP5) Homo sapiens Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which may mediate cytoskeleton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],