

RHG30 Polyclonal Antibody

Cat No: HR1AP11274

For research use only

Overview

Product Name	RHG30 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	121kDa
GeneID?Human?	ARHGAP30
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.,similarity:Contains 1 Rho-GAP domain.,