

RHG28 Polyclonal Antibody

Cat No: HR1AP11273

For research use only

Overview

Product Name	RHG28 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	ARHGAP28 KIAA1314
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 Rho-GAP domain.,tissue specificity:Expressed in testis. Expressed at moderate level in kidney and ovary, and weakly expressed in spleen and skeletal muscle.,