

RHG25 Polyclonal Antibody

Cat No: HR1AP11272

For research use only

Overview

Product Name	RHG25 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	ARHGAP25 KIAA0053
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Rho GTPase activating protein 25(ARHGAP25) Homo sapiens ARHGAPs, such as ARHGAP25, encode negative regulators of Rho GTPases (see ARHA; MIM 165390), which are implicated in actin remodeling, cell polarity, and cell migration (Katoh and Katoh, 2004 [PubMed 15254788]).[supplied by OMIM, Mar 2008],