

RHG24 Polyclonal Antibody

Cat No: HR1AP11271

For research use only

Overview

Product Name	RHG24 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	82kDa
GeneID?Human?	ARHGAP24 FILGAP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Rho GTPase activating protein 24(ARHGAP24) Homo sapiens This gene encodes a Rho-GTPase activating protein, which is specific for the small GTPase family member Rac. Binding of the encoded protein by filamin A targets it to sites of membrane protrusion, where it antagonizes Rac. This results in suppression of lamellae formation and promotion of retraction to regulate cell polarity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016].