

RHG20 Polyclonal Antibody

Cat No: HR1AP11270

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | RHG20 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 131kDa |
| GeneID?Human? | ARHGAP20 KIAA1391 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | Rho GTPase activating protein 20(ARHGAP20) Homo sapiens The protein encoded by this gene is an activator of RHO-type GTPases, transducing a signal from RAP1 to RHO and impacting neurite outgrowth. [provided by RefSeq, Sep 2016], |