

RBBP6 Polyclonal Antibody

Cat No: HR1AP11267

For research use only

Overview

Product Name	RBBP6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	197kDa
GeneID?Human?	RBBP6 P2PR PACT RBQ1 My038
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RB binding protein 6, ubiquitin ligase(RBBP6) Homo sapiens The retinoblastoma tumor suppressor (pRB) protein binds with many other proteins. In various human cancers, pRB suppresses cellular proliferation and is inactivated. Cell cycle-dependent phosphorylation regulates the activity of pRB. This gene encodes a protein which binds to underphosphorylated but not phosphorylated pRB. Multiple alternatively spliced transcript variants that encode different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],</p>