



RFPL2 Polyclonal Antibody

Cat No: HR1AP11265

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | RFPL2 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 41kDa |
| GeneID?Human? | RFPL2 RNF79 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 RING-type zinc finger.,tissue specificity:Seems to be expressed in prostate and less abundantly in adult brain, fetal liver, and fetal kidney., |