

# RSF1 Polyclonal Antibody

Cat No: HR1AP11263

For research use only

## Overview

Product Name	RSF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	158kDa
GeneID?Human?	RSF1 HBXAP XAP8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	remodeling and spacing factor 1(RSF1) Homo sapiens This gene encodes a nuclear protein that interacts with hepatitis B virus X protein (HBX) and facilitates transcription of hepatitis B virus genes by the HBX transcription activator, suggesting a role for this interaction in the virus life cycle. This protein also interacts with SNF2H protein to form the RSF chromatin-remodeling complex, where the SNF2H subunit functions as the nucleosome-dependent ATPase, and this protein as the histone chaperone. [provided by RefSeq, Sep 2011],