

# RGS22 Polyclonal Antibody

Cat No: HR1AP11261

For research use only

## Overview

Product Name	RGS22 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	139kDa
GeneID?Human?	RGS22
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits thereby driving them into their inactive GDP-bound form.,similarity:Contains 2 RGS domains.,