

RRAS2 Polyclonal Antibody

Cat No: HR1AP11249

For research use only

Overview

Product Name	RRAS2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RRAS2 TC21
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>related RAS viral (r-ras) oncogene homolog 2(RRAS2) Homo sapiens This gene encodes a member of the R-Ras subfamily of Ras-like small GTPases. The encoded protein associates with the plasma membrane and may function as a signal transducer. This protein may play an important role in activating signal transduction pathways that control cell proliferation. Mutations in this gene are associated with the growth of certain tumors. Pseudogenes of this gene are found on chromosomes 1 and 2. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Apr 2010],</p>