

RAB9B Polyclonal Antibody

Cat No: HR1AP11247

For research use only

Overview

Product Name	RAB9B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RAB9B RAB9L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB9B, member RAS oncogene family(RAB9B) Homo sapiens This gene encodes a member of a subfamily of RAS small guanosine triphosphate (GTP)-binding proteins that regulate membrane trafficking. The encoded protein may be involved in endosome-to-Golgi transport. [provided by RefSeq, Jan 2010],