

RAB9A Polyclonal Antibody

Cat No: HR1AP11246

For research use only

Overview

Product Name	RAB9A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RAB9A RAB9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Involved in the transport of proteins between the endosomes and the trans Golgi network.,similarity:Belongs to the small GTPase superfamily. Rab family.,