

RALB Polyclonal Antibody

Cat No: HR1AP11248

For research use only

Overview

Product Name	RALB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RALB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAS like proto-oncogene B(RALB) Homo sapiens This gene encodes a GTP-binding protein that belongs to the small GTPase superfamily and Ras family of proteins. GTP-binding proteins mediate the transmembrane signaling initiated by the occupancy of certain cell surface receptors. [provided by RefSeq, Jul 2008],