

# RAB6B Polyclonal Antibody

Cat No: HR1AP11243

For research use only

## Overview

Product Name	RAB6B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RAB6B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Seems to have a role in retrograde membrane traffic at the level of the Golgi complex.,similarity:Belongs to the small GTPase superfamily. Rab family.,subunit:Interacts with RAB6KIFL.,tissue specificity:Predominantly expressed in brain.,