

RB40L Polyclonal Antibody

Cat No: HR1AP11239

For research use only

Overview

Product Name	RB40L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	RAB40AL RLGP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB40A, member RAS oncogene family-like(RAB40AL) Homo sapiens This gene encodes a member of the Rab40 subfamily of Rab small GTP-binding proteins that contains a C-terminal suppressors of cytokine signaling box. Disruptions in this gene are associated with Duchenne muscular dystrophy. [provided by RefSeq, Apr 2010],