

RAB43 Polyclonal Antibody

Cat No: HR1AP11241

For research use only

Overview

Product Name	RAB43 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RAB43 RAB41
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Widely expressed in brain, testis, lung, heart, ovary, colon, kidney, uterus and spleen but not in liver.,