

RB40A Polyclonal Antibody

Cat No: HR1AP11238

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | RB40A Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 30kDa |
| GeneID?Human? | RAB40A |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | RAB40A, member RAS oncogene family(RAB40A) Homo sapiens This gene encodes a member of the Rab40 subfamily of Rab small GTP-binding proteins that contains a C-terminal suppressors of cytokine signaling box. [provided by RefSeq, Apr 2010], |