

## RB39B Polyclonal Antibody

Cat No: HR1AP11237

For research use only

### Overview

Product Name	RB39B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RAB39B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB39B, member RAS oncogene family(RAB39B) Homo sapiens This gene encodes a member of the Rab family of proteins. Rab proteins are small GTPases that are involved in vesicular trafficking. Mutations in this gene are associated with X-linked mental retardation. [provided by RefSeq, Aug 2013],