

RB39A Polyclonal Antibody

Cat No: HR1AP11236

For research use only

Overview

Product Name	RB39A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RAB39A RAB39
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in vesicular trafficking.,similarity:Belongs to the small GTPase superfamily. Rab family.,