

RAB28 Polyclonal Antibody

Cat No: HR1AP11231

For research use only

Overview

Product Name	RAB28 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	RAB28
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB28, member RAS oncogene family(RAB28) Homo sapiens This gene encodes a member of the Rab subfamily of Ras-related small GTPases. The encoded protein may be involved in regulating intracellular trafficking. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 9 and X. [provided by RefSeq, Apr 2009],