

RAB32 Polyclonal Antibody

Cat No: HR1AP11233

For research use only

Overview

Product Name	RAB32 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	RAB32
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB32, member RAS oncogene family(RAB32) Homo sapiens The protein encoded by this gene anchors the type II regulatory subunit of protein kinase A to the mitochondrion and aids in mitochondrial fission. The encoded protein also appears to be involved in autophagy and melanosome secretion. Variations in this gene may be linked to leprosy. [provided by RefSeq, Dec 2015],