

## RB27A Polyclonal Antibody

Cat No: HR1AP11230

For research use only

### Overview

Product Name	RB27A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	RAB27A RAB27
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB27A, member RAS oncogene family(RAB27A) Homo sapiens The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. The protein is membrane-bound and may be involved in protein transport and small GTPase mediated signal transduction. Mutations in this gene are associated with Griscelli syndrome type 2. Alternative splicing occurs at this locus and four transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008],