

RB22A Polyclonal Antibody

Cat No: HR1AP11228

For research use only

Overview

Product Name	RB22A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	RAB22A RAB22
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB22A, member RAS oncogene family(RAB22A) Homo sapiens The protein encoded by this gene is a member of the RAB family of small GTPases. The GTP-bound form of the encoded protein has been shown to interact with early-endosomal antigen 1, and may be involved in the trafficking of and interaction between endosomal compartments. [provided by RefSeq, Jul 2008],