



RAB17 Polyclonal Antibody

Cat No: HR1AP11224

For research use only

Overview

Product Name	RAB17 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RAB17
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB17, member RAS oncogene family(RAB17) Homo sapiens The Rab subfamily of small GTPases plays an important role in the regulation of membrane trafficking. RAB17 is an epithelial cell-specific GTPase (Lutcke et al., 1993 [PubMed 8486736]).[supplied by OMIM, Oct 2009],