

RAB1A Polyclonal Antibody

Cat No: HR1AP11226

For research use only

Overview

Product Name	RAB1A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RAB1A RAB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RAB1A, member RAS oncogene family(RAB1A) Homo sapiens This gene encodes a member of the Ras superfamily of GTPases. Members of the gene family cycle between inactive GDP-bound and active GTP-bound forms. This small GTPase controls vesicle traffic from the endoplasmic reticulum to the Golgi apparatus. Multiple alternatively spliced transcript variants have been identified for this gene which encode different protein isoforms. [provided by RefSeq, Oct 2008],</p>