

RAB19 Polyclonal Antibody

Cat No: HR1AP11225

For research use only

Overview

Product Name	RAB19 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RAB19 RAB19B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the small GTPase superfamily. Rab family.,