

RAB15 Polyclonal Antibody

Cat No: HR1AP11223

For research use only

Overview

Product Name	RAB15 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RAB15
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May act in concert with RAB3A in regulating aspects of synaptic vesicle membrane flow within the nerve terminal.,similarity:Belongs to the small GTPase superfamily. Rab family.,subunit:The GTP bound form of RAB15 interacts with REP15.,