

# RRAGD Polyclonal Antibody

Cat No: HR1AP11221

For research use only

## Overview

Product Name	RRAGD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	RRAGD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Ras related GTP binding D(RRAGD) Homo sapiens RRAGD is a monomeric guanine nucleotide-binding protein, or G protein. By binding GTP or GDP, small G proteins act as molecular switches in numerous cell processes and signaling pathways.[supplied by OMIM, Apr 2004],