

## RAC3 Polyclonal Antibody

Cat No: HR1AP11220

For research use only

### Overview

Product Name	RAC3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	RAC3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ras-related C3 botulinum toxin substrate 3 (rho family, small GTP binding protein Rac3)(RAC3) Homo sapiens The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2015],</p>