

RGF1A Polyclonal Antibody

Cat No: HR1AP11215

For research use only

Overview

Product Name	RGF1A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	RASGEF1A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Guanine nucleotide exchange factor (GEF) for KRAS, HRAS, and NRAS (in vitro). Plays a role in cell migration.,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-GEF domain.,tissue specificity:Detected in brain and spinal cord. Highly expressed in a number of intrahepatic cholangiocarcinoma tissue biopsies.,