

RGF1B Polyclonal Antibody

Cat No: HR1AP11216

For research use only

Overview

Product Name	RGF1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	RASGEF1B GPIG4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Guanine nucleotide exchange factor (GEF) for Ras family proteins (in vitro).,induction:Up-regulated in monocytes stimulated with T.cruzi GPI-anchored mucins or bacterial lipopolysaccharide.,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-GEF domain.,