

RAPH1 Polyclonal Antibody

Cat No: HR1AP11213

For research use only

Overview

Product Name	RAPH1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	137kDa
GeneID?Human?	RAPH1 ALS2CR18 ALS2CR9 KIAA1681 LPD PREL2 RMO1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>Ras association (RalGDS/AF-6) and pleckstrin homology domains 1(RAPH1) Homo sapiens This gene encodes a protein that belongs to the Mig10/Rap1-interacting adaptor molecule/Lamellipodin family of adapter proteins, which function in cell migration. Members of this family contain pleckstrin-homology domains, Ras-association domains, and proline-rich C-termini. The protein encoded by this gene regulates actin dynamics through interaction with Ena/Vasodilator proteins as well as direct binding to filamentous actin to regulate actin network assembly. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2016],</p>