

IQGA2 Polyclonal Antibody

Cat No: HR1AP11211

For research use only

Overview

Product Name	IQGA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	173kDa
GeneID?Human?	IQGAP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>IQ motif containing GTPase activating protein 2(IQGAP2) Homo sapiens This gene encodes a member of the IQGAP family. The protein contains three IQ domains, one calponin homology domain, one Ras-GAP domain and one WW domain. It interacts with components of the cytoskeleton, with cell adhesion molecules, and with several signaling molecules to regulate cell morphology and motility. [provided by RefSeq, Jul 2008].</p>