



## GDS1 Polyclonal Antibody

Cat No: HR1AP11207

For research use only

### Overview

|                       |  |
|-----------------------|--|
| Product Name          | GDS1 Polyclonal Antibody   |
| Source                | Rabbit   |
| Applications          | WB,ELISA   |
| Species Reactivity    | Human  |
| Recommended Dilutions |  |
| Immunogen             |  |
| Species               | Rabbit   |
| Storage               | -20°C/1 year   |
| Isotype               |  |
| Clonality             |  |
| Concentration         | 1 mg/ml  |
| Observed band         | 66kDa  |
| GeneID?Human?         | RAP1GDS1   |
| Human Swiss-Prot No.  |  |
| Cellular localization |  |
| Alternative Names     |  |
| Background            | Rap1 GTPase-GDP dissociation stimulator 1(RAP1GDS1) Homo sapiens The smg GDP dissociation stimulator (smgGDS) protein is a stimulatory GDP/GTP exchange protein with GTPase activity (Riess et al., 1993 [PubMed 8262526]).[supplied by OMIM, Feb 2010], |