

REPS2 Polyclonal Antibody

Cat No: HR1AP11201

For research use only

Overview

Product Name	REPS2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	72kDa
GeneID?Human?	REPS2 POB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RALBP1 associated Eps domain containing 2(REPS2) Homo sapiens The product of this gene is part of a protein complex that regulates the endocytosis of growth factor receptors. The encoded protein directly interacts with a GTPase activating protein that functions downstream of the small G protein Ral. Its expression can negatively affect receptor internalization and inhibit growth factor signaling. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],</p>