



GNDS Polyclonal Antibody

Cat No: HR1AP11198

For research use only

Overview

Product Name	GNDS Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	RALGDS KIAA1308 RGF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ral guanine nucleotide dissociation stimulator(RALGDS) Homo sapiens Guanine nucleotide dissociation stimulators (GDSs, or exchange factors), such as RALGDS, are effectors of Ras-related GTPases (see MIM 190020) that participate in signaling for a variety of cellular processes.[supplied by OMIM, Nov 2010],