

# RSPH1 Polyclonal Antibody

Cat No: HR1AP11196

For research use only

## Overview

Product Name	RSPH1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	RSPH1 TSA2 TSGA2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>radial spoke head 1 homolog(RSPH1) Homo sapiens This gene encodes a male meiotic metaphase chromosome-associated acidic protein. This gene is expressed in tissues with motile cilia or flagella, including the trachea, lungs, airway brushings, and testes. Mutations in this gene result in primary ciliary dyskinesia-24. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2014],</p>