

# RGL1 Polyclonal Antibody

Cat No: HR1AP11199

For research use only

## Overview

Product Name	RGL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GeneID?Human?	RGL1 KIAA0959 RGL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Probable guanine nucleotide exchange factor.,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-associating domain.,similarity:Contains 1 Ras-GEF domain.,subunit:Interacts with Ras.,tissue specificity:Expressed in a wide variety of tissues with strong expression being seen in the heart, brain, kidney, spleen and testis.,