

## RBL2B Polyclonal Antibody

Cat No: HR1AP11193

For research use only

### Overview

Product Name	RBL2B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	RABL2B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB, member of RAS oncogene family-like 2B(RABL2B) Homo sapiens The RABL2B protein is a member of the RAB gene family which belongs to the RAS GTPase superfamily. RABL2B is located within a subtelomeric region of 22q13.3. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2008],