

DEN5A Polyclonal Antibody

Cat No: HR1AP11191

For research use only

Overview

Product Name	DEN5A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	141kDa
GeneID?Human?	DENND5A KIAA1091 RAB6IP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in Rab6-mediated GTPase signaling.,similarity:Belongs to the RAB6IP1 family.,similarity:Contains 1 dDENN domain.,similarity:Contains 1 DENN domain.,similarity:Contains 1 PLAT domain.,similarity:Contains 1 uDENN domain.,similarity:Contains 2 RUN domains.,subunit:Probably interacts with RAB6A bound to GTP.,