

RBGP1 Polyclonal Antibody

Cat No: HR1AP11189

For research use only

Overview

Product Name	RBGP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	117kDa
GeneID?Human?	RABGAP1 HSPC094
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:May act as a GTPase-activating protein of RAB6A. May play a role in microtubule nucleation by centrosome. May participate in a RAB6A-mediated pathway involved in the metaphase-anaphase transition.,sequence caution:Contaminating sequence. Potential poly-A sequence.,sequence caution:Intron retention.,similarity:Contains 1 PID domain.,similarity:Contains 1 Rab-GAP TBC domain.,subcellular location:Predominantly cytosolic but also associated with the centrosome.,subunit:Interacts with RAB6A and tubulin gamma.,</p>