

NSUN5 Polyclonal Antibody

Cat No: HR1AP11184

For research use only

Overview

Product Name	NSUN5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	NSUN5 WBSCR20 WBSCR20A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>NOP2/Sun RNA methyltransferase family member 5(NSUN5) Homo sapiens This gene encodes a member of an evolutionarily conserved family of proteins that may function as methyltransferases. This gene is located in a larger region of chromosome 7 that is deleted in Williams-Beuren syndrome, a multisystem developmental disorder. There are two pseudogenes for this gene located in the same region of chromosome 7. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2013].</p>