

ZPI Polyclonal Antibody

Cat No: HR1AP11180

For research use only

Overview

Product Name	ZPI Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	SERPINA10 ZPI UNQ707/PRO1358
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	serpin family A member 10(SERPINA10) Homo sapiens The protein encoded by this gene belongs to the serpin family. It is predominantly expressed in the liver and secreted in plasma. It inhibits the activity of coagulation factors Xa and XIa in the presence of protein Z, calcium and phospholipid. Mutations in this gene are associated with venous thrombosis. Alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, May 2010].