

PCDH7 Polyclonal Antibody

Cat No: HR1AP11182

For research use only

Overview

Product Name	PCDH7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	117kDa
GeneID?Human?	PCDH7 BHPCDH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>protocadherin 7(PCDH7) Homo sapiens This gene belongs to the protocadherin gene family, a subfamily of the cadherin superfamily. The gene encodes a protein with an extracellular domain containing 7 cadherin repeats. The gene product is an integral membrane protein that is thought to function in cell-cell recognition and adhesion. Alternative splicing yields isoforms with unique cytoplasmic tails. [provided by RefSeq, Jul 2008],</p>