

S10AE Polyclonal Antibody

Cat No: HR1AP11174

For research use only

Overview

Product Name	S10AE Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	S100A14 S100A15
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	S100 calcium binding protein A14(S100A14) Homo sapiens This gene encodes a member of the S100 protein family which contains an EF-hand motif and binds calcium. The gene is located in a cluster of S100 genes on chromosome 1. Levels of the encoded protein have been found to be lower in cancerous tissue and associated with metastasis suggesting a tumor suppressor function (PMID: 19956863, 19351828). [provided by RefSeq, Dec 2011],