

PPME1 Polyclonal Antibody

Cat No: HR1AP11172

For research use only

Overview

Product Name	PPME1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	PPME1 PME1 PP2593 PRO0750
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>protein phosphatase methylesterase 1(PPME1) Homo sapiens This gene encodes a protein phosphatase methylesterase localized to the nucleus. The encoded protein acts on the protein phosphatase-2A catalytic subunit and supports the ERK pathway through dephosphorylation of regulatory proteins. It plays a role in malignant glioma progression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2012],</p>