

PPM1H Polyclonal Antibody

Cat No: HR1AP11169

For research use only

Overview

Product Name	PPM1H Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	PPM1H ARHCL1 KIAA1157 URCC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:A phosphoprotein + H(2)O = a protein + phosphate.,similarity:Belongs to the PP2C family.,similarity:Contains 1 PP2C-like domain.,